

**SnaB I**

Cat. No.	용량	농도
DYR1950	500 units	5 units/μl
DYR1952	1,000 units	5 units/μl
DYR1954	2,500 units	5 units/μl
DYR1956	2,500 units	25 units/μl

**◆ 제품구성**

- SnaB I
- 10X DY Buffer IV
- 10X FastCut Buffer
- Sterile water
- Dyne 6X DNA Loading Buffer ver.2

**◆ Source**

- Sphaerotilus natans*

**◆ Quality control**

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
- Extreme purity assay

**◆ 인식부위**



Single letter code	
W = A or T	S = C or V = A or C or G
N = A or C or G or TG	M = A or C
K = G or T	R = A or G
Y = C or T	B = C or G or T
D = A or G or T	H = A or C or T

**◆ 보관온도**

- 20°C

**◆ Heat inactivation**

- 80°C, 20분.

**◆ Unit정의**

- 1 unit은 박테리오파지 DNA(Hind III digestion) 1 μg을 50 μl 반응물로 37°C에서 1시간 동안 완전히 분해하는데 필요한 효소의 양이다.

**◆ Buffer별 상대적 활성도**

I	II	III	IV	FastCut
100%	75%	25%	100%	100%

**◆ Methylation effect**

Methylation	dam	dcm	CpG
Cleavage	Cleavage	Cleavage	No Cleavage

**◆ 주의사항**

- CpG 메틸화(methylation)은 mammalian genomicDNA 절단을 저해한다. 반응은 3시간을 넘지 않도록 한다. 과량의 효소, 고농도 glycerol(>5%) 또는 높은 pH(>8.0)의 배양상태에서는 비특이적인 활성이 나타난다.

**◆ 표준반응 조건**

- Normal Protocol

Component	농도	Volume
Substrate DNA	1 μg	X μl
10X DY Buffer IV	1 X	5 μl
SnaB I		Substrate dependent
Sterile water		Up to 50 μl

\* Incubate at 37°C for 1 hr

- Fast Protocol

Component	농도	Volume
Substrate DNA	1 μg	X μl
10X FastCut Buffer	1 X	5 μl
SnaB I	5 unit	1 μl
Sterile water		Up to 50 μl

\* Incubate at 37°C for 15 min